

# SimDis Expert

## Version 10

The Evolution of Excellence Continues



# The 'Gold' Standard

**SimDis® Expert** is considered the gold standard for GC simulated distillation software. Its rich set of capabilities and very intuitive graphical user interface have made the 'GC data to decision' workflow simple and straightforward since its introduction 20 years ago. And, with Version 10, SimDis® Expert is more powerful and flexible than ever.

SimDis® Expert meets or exceeds all of the requirements of existing and proposed standard simulated distillation test methods. It's fully compatible with all of today's chromatography data systems and can be tailored to your specific needs including support for selective detectors including FPD, PFPD, SCD, NCD and even MS. It's available for installation on a PC as a client or on a network server.

Like all of our software, SimDis® Expert is fully backed by a highly responsive support team with an unrivaled level of expertise in all of the GC and GC-MS based test methods used in the hydrocarbon processing industry today.

## Key Capabilities of SimDis® Expert 10

- Supports all current and proposed standard test methods
- Utilizes a graphical user interface with straightforward analysis workflow
- Simulated Wax Analysis Report [Similar to ASTM D5442] **(NEW)**
- Cut points can be calculated as average molecular weight **(NEW)**
- Configurable automatic baseline selection [Before, After or Both] **(NEW)**
- The test method configuration screen now includes detailed information about the scope of the selected method, GC conditions, recommended standards, reference materials and consumables **(NEW)**
- Reports can be easily customized; ex. place boiling point curve and cut point table on a single page **(NEW)**
- Auto-configuring 'universal' software driver to support all GC chromatography data systems **(NEW)**
- Supports selective detectors (SCD, FPD, PFPD, AED) (Advanced Version) **(NEW)**
- Optional MS module for aromatics and saturates analysis (Advanced Version) **(NEW)**
- D86 and D1160 correlations
- Crude oil blending simulation model including aromatics and saturates analysis (ASA) (Advanced Version Only)
- Built in Merge Expert™ for DHA 'front end' correction of crude oil samples (Advanced Version)
- Full range of standard calculations: MOV, NOACK, Reid Vapor Pressure (RVP), CETANE index, sample average molecular weight and others.
- Built in peak integrator with expanded integration parameters to simplify data reprocessing
- Full range of standard report types as well as the ability to create your own custom reports
- Fully customizable cut points
- Peak skew and column resolution determination based on multiple peaks
- User access levels can be assigned to meet data security requirements
- Results can be automatically exported as delimited ASCII/CSV (LIMS), to Microsoft Excel, OpenOffice.
- Screen views can be exported to the Windows 'clipboard' or as Windows Enhanced Metafiles
- 'Intelligent' adjustment of the peak/elution detection algorithm; it 'learns' based on user defined elution marks
- Custom calculation formula editor for 'Interactive' distillation; ex. display the percentage distilled and temperature at any point of the chromatogram
- Special software programs can be executed automatically in the event of a QC failure (ex. send text message) **and much more.....**



# Supports All Test Methods

## SimDis® Expert is Always Up to Date

SimDis® Expert is routinely updated to ensure it continually meets or exceeds our client's simulated distillation analysis requirements. Standard test methods are added as they're approved or once they've been deemed a proposed test method by the standards testing organization committee.

With the release of Version 10, GC operational conditions for each test method and a list of calibration standards and reference materials appropriate for the selected method is shown on the software's method setup window'. (see below)

**Available Test Methods**


**GC Column and Conditions**

**Recommended Standards & Reference Materials**

Method	Column	Carrier	Carrier Flow	Initial column temperature C	Final column temperature C	Oven rate	Final hold time	Detector	Detector temperature	Injector	Injector temperature	Sample size	Sample concentration
Method A	7.5 m x 0.53 mm x 1.5 µm (DB-1)	Nitrogen	30 mL/min	40	340	10 °C/min	0	FID	350 °C	PTV	340 °C	0.5 µL	25 (mass %)
Method-1B	5 m x 0.53 mm x 0.88 µm (HP-1)	Helium	12 mL/min	35	350	10 °C/min	0	FID	380 °C	COC	≈ 3 °C > oven	1 µL	2 (mass %)
Method-1C													

PartNo	Description
SD-060-05	Retention time standard without Ethanol
SD-060	Retention time standard with Ethanol
SD-007-RF	Response factor standard
SD-008-RT	Retention time standard

# Basic Capabilities




An example of a retention time calibration chromatogram with associated boiling point curve for ASTM D2887 (above) and the corresponding boiling point table with automatic skewness fit. (right)

Component	Time	BP(F)	Skewness
n-C5	0.198	97.0	
n-C6	0.332	156.0	
n-C7	0.613	209.2	
n-C8	1.126	259.0	
n-C9	1.872	304.0	
n-C10	2.763	345.4	
n-C11	3.696	385.0	
n-C12	4.618	421.3	
n-C13	5.495	455.7	
n-C14	6.337	489.0	1.07
n-C15	7.130	520.0	
n-C16	7.886	549.0	
n-C17	8.602	576.0	
n-C18	9.285	601.0	
n-C19	9.937	626.2	
n-C20	10.558	651.0	
n-C24	12.795	736.0	
n-C28	14.720	808.3	
n-C32	16.400	870.3	
n-C36	17.897	924.8	
n-C40	19.237	971.6	
n-C44	20.447	1013.0	

File	D:\Work\sde\Demo Files\Regular\D2887_Calib.SSF
Sample Id	SS3E025T
LIMS Id	
Parameter File	D2887
Injection Date	1/5/2001 5:13:50 PM
Operator	Demo Chemist

# Advanced Capabilities




## Powerful Peak Integration

With the built-in integrator, making adjustments to peak integration is simple. There is no need to switch back and forth between the GC data acquisition system used to acquire the signals and SimDis Expert. You can stay within the SimDis Expert environment for all your reprocessing needs. This is especially useful handling retention time calibration standards containing heavy carbons such as n-C<sub>110</sub> and n-C<sub>120</sub>

## Solvent Quenching Tool

The software has an algorithm to automatically detect the quenching window which will change as the column ages. As such, the beginning and end of the solvent peak can be easily adjusted using the interactive graphical interface.



## Automatic Baseline Determination


When loading a signal for interactive analysis, the software can automatically determine the most suitable baseline.

When the software is set up to operate in an unattended manner (automatic execution after each sample injection) the user can define whether samples should be corrected using the baseline obtained prior or using one acquired at a later time.

The screenshot shows a data table with the following columns: Signal, Sample ID, Parameter, Acq. Method, Mode, Channel, Blank, and QC. A warning message is displayed in the top right corner: "Warning: Sample or solvent amount cannot be zero. This method needs to calculate a recovery." Below the table, there are options for automatic baseline selection, including "None", "Before sample", "After sample", and "Before or after sample".

Signal	Sample ID	Parameter	Acq. Method	Mode	Channel	Blank	QC
D:\Work\pde\Demo Files\RegAir\ASTM7169_RESID_310A.DSP	SI006_Y1	ASTM7169_Y10	D7169_RESIDUE	Standard	PD(S)		1
D:\Work\pde\Demo Files\RegAir\ASTM7169_RESID_310A.DSP	RESID5_V4L38A	ASTM7169_Y10	D7169_RESIDUE	Sample	PD(S)		1

## More Advanced Capabilities



### Selective Detector Capability

Seamlessly process signals from selective detector(s) in addition to the FID or in place of it. Extend traditional simulated distillation to sulfur (ex. ASTM D7807), nitrogen and other elements through the use of selective detectors such as SCD, NCD, FPD, PFPD, AED and MSD.

### Advanced and Custom Calculations


Extend your analysis results through the use of additional calculations and/or formulas defined by the user.

The 'Additional Calculations' dialog box is shown with the following settings:

- Correlation:**  Automatic,  D-86,  D-1160,  None
- Calculated Flash Point (ASTM D-7215):**
  - CFP-D56
  - CFP-D93
  - CFP-D3828
  - Maximum MSPE<sub>x</sub> allowed:  [F]
- Other:**
  - Reid Vapor Pressure
  - Vapor / Liquid Ratio
  - Volume Percent Conversion
    - Polynomial model (Requires crude's density)
    - ASTM X:XXXX
  - Volatility at  [F]
  - Cetane Index by ASTM D-976
  - Cetane Index by ASTM D-4737
  - NOACK Volatility (DIN 51581-2) @  [F]

Buttons: Ok, Help, Cancel

# Crude Oil Blending Simulator



## Available Blending Simulation Model for Crude Oil and Fully Eluting Materials

Save time and money by simulating a gravimetric blend comprised of up to 20 different samples. With a single click of the mouse the software creates a virtual signal and properties of the simulated gravimetric blend.

Accuracy of the simulated blend (green box) vs. the actual blend (purple box) below

Cut(F)	%Off	%Off	(-)
( SET, 350.0 )	8.08	8.19	-0.11
( 350.0, 400.0 )	2.35	2.41	-0.07
( 400.0, 450.0 )	2.88	2.99	-0.10
( 450.0, 500.0 )	3.52	3.59	-0.06
( 500.0, 550.0 )	4.00	3.93	0.07
( 550.0, 600.0 )	4.77	4.85	-0.08
( 600.0, 650.0 )	4.69	4.78	-0.08
( 650.0, 700.0 )	4.66	4.58	0.08
( 700.0, 750.0 )	4.59	4.57	0.02
( 750.0, 800.0 )	5.19	5.18	0.01
( 800.0, 850.0 )	5.42	5.31	0.11
( 850.0, 900.0 )	4.80	4.73	0.07
( 900.0, 950.0 )	4.45	4.42	0.04
( 950.0, 1000.0 )	4.17	4.15	0.02
( 1000.0, 1050.0 )	4.05	4.04	0.01
( 1050.0, 1100.0 )	4.54	4.55	-0.01
( 1100.0, 1150.0 )	5.09	5.10	-0.01
( 1150.0, 1200.0 )	4.85	4.79	0.06
( 1200.0, 1250.0 )	3.95	3.83	0.12
( 1250.0, 1300.0 )	3.58	3.43	0.15
( 1300.0, EET )	3.09	2.87	0.22

Sample Id	Start Time	End Time	%Recovery
100Hamaca\$A	0.000	29.700	92.72
50/50HamPtbra\$A	0.000	29.700	92.28

%Off	BP(F)	BP(F)	(-)
68.00	1054.5	1053.5	1.0
69.00	1065.5	1064.5	1.0
70.00	1076.8	1075.8	1.0
71.00	1087.8	1086.8	1.0
72.00	1098.4	1097.4	1.0
73.00	1108.5	1107.5	1.0
74.00	1118.6	1117.5	1.0
75.00	1128.5	1127.5	1.0
76.00	1138.1	1137.1	1.0
77.00	1147.6	1146.6	1.1
78.00	1157.2	1156.2	1.1
79.00	1167.1	1166.1	1.0
80.00	1177.0	1176.1	0.9
81.00	1187.7	1186.8	0.8
82.00	1198.9	1198.3	0.6
83.00	1209.9	1209.6	0.3
84.00	1221.7	1221.7	0.0
85.00	1235.9	1236.4	-0.5
86.00	1249.2	1250.3	-1.0
87.00	1264.2	1265.5	-1.4
88.00	1278.0	1279.6	-1.6
89.00	1291.8	1294.2	-2.3
90.00	1304.8	1308.3	-3.5
91.00	1319.5	1323.9	-4.4
92.00	1334.3	1340.5	-6.1



  

Sample Id	Start Time	End Time	%Recovery
100Hamaca\$A	0.000	29.700	92.72
50/50HamPtbra\$A	0.000	29.700	92.28

# Powerful Yet Exceptionally Easy to Use

## Exceptionally Easy to Use


Highly intuitive graphical user interface makes it very easy to review data, compare signals, generate plots and apply custom calculations.



## Boiling Point Curve 'Front End Correction' Using DHA Data

The boiling point curve for crude oil produced by high temperature simulated distillation can be corrected through the use of the included Merge Expert™ software (Advanced Version) to process boiling point data produced by a GC running any of the standard detailed hydrocarbon analysis (DHA) test methods

%Off	BP(F)	BP(C)	(-)
IBP	93.2	-23.2	116.4
1.00	94.8	40.7	54.1
2.00	101.3	83.5	17.8
3.00	151.8	127.3	24.5
4.00	165.3	150.5	14.8
5.00	190.9	181.5	9.4
6.00	211.6	199.7	11.9
7.00	236.7	223.7	13.0
8.00	264.4	248.0	16.4
9.00	287.5	275.2	12.2
10.00	314.9	299.5	15.4
11.00	336.4	336.1	0.3
12.00	359.9	362.9	-2.9
13.00	384.3	384.3	0.0
14.00	405.0	405.0	0.0
15.00	423.2	423.2	0.0
16.00	439.6	439.6	0.0
17.00	456.0	456.0	0.0
18.00	471.8	471.8	0.0
19.00	486.1	486.1	0.0
20.00	499.8	499.8	0.0
21.00	512.3	512.3	0.0
22.00	525.7	525.7	0.0
23.00	540.2	540.2	0.0
24.00	553.4	553.4	0.0
25.00	565.4	565.4	0.0
26.00	577.1	577.1	0.0





The one source for all your simulated distillation needs

[www.separationsystems.com](http://www.separationsystems.com)

## About Us

We are a minority owned business located in Gulf Breeze, Florida. We offer GC and GC-MS based analysis systems, application software, consumables, support and training for petroleum refining, bio-fuels and petrochemical applications. Our systems are comprehensive in nature and include a GC or GC-MS, our own specialized hardware and software, reference & calibration standards, consumables, training and support.

While the majority of our systems are designed to meet the international standard testing method requirements (ex. ASTM, EN, ISO), we also design systems for special requirements including custom software.

## Headquarters

Separation System, Inc.  
100 Nightingale Lane  
Gulf Breeze, Florida 32651  
USA

Telephone: 1-800-340-3302  
Email: [sales@separationsystems.com](mailto:sales@separationsystems.com)  
Web: [www.separationsystems.com](http://www.separationsystems.com)